

file caplus
COST IN U.S. DOLLARS

SINCE FILE	TOTAL
ENTRY	SESSION
48.17	69.02

FULL ESTIMATED COST

FILE 'CAPLUS' ENTERED AT 16:40:26 ON 24 DEC 2003
USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.
COPYRIGHT (C) 2003 AMERICAN CHEMICAL SOCIETY (ACS)

Copyright of the articles to which records in this database refer is held by the publishers listed in the PUBLISHER (PB) field (available for records published or updated in Chemical Abstracts after December 26, 1996), unless otherwise indicated in the original publications. The CA Lexicon is the copyrighted intellectual property of the American Chemical Society and is provided to assist you in searching databases on STN. Any dissemination, distribution, copying, or storing of this information, without the prior written consent of CAS, is strictly prohibited.

FILE COVERS 1907 - 24 Dec 2003 VOL 139 ISS 26
FILE LAST UPDATED: 23 Dec 2003 (20031223/ED)

This file contains CAS Registry Numbers for easy and accurate substance identification.

=> s (paral?(w) mass (w)spectr?)

321747 PARAL?
788655 MASS
71248 MASSES
826335 MASS
(MASS OR MASSES)

2285432 SPECTR?

L10 16 (PARAL?(W) MASS (W)SPECTR?)

=> d bib,abs 1-16

L10 ANSWER 1 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN

AN 2003:289440 CAPLUS

DN 139:193345

TI A generic assay for phosphate-consuming or -releasing enzymes coupled on-line to liquid chromatography for lead finding in natural products
AU Schenk, T.; Appels, N. M. G. M.; van Elswijk, D. A.; Irth, H.; Tjaden, U. R.; van der Greef, J.

CS Kiadis B.V., Leiden, 2333 CA, Neth.

SO Analytical Biochemistry (2003), 316(1), 118-126

CODEN: ANBCA2; ISSN: 0003-2697

PB Elsevier Science

DT Journal

LA English

AB A generic continuous-flow assay for phosphate-consuming or -releasing enzymes coupled online to liq. chromatog. (LC) has been developed. Operating the LC-biochem. assay in combination with mass spectrometry allows the fast detection and identification of inhibitors of these enzymes in complex mixts. The assay is based on the detection of phosphate, released by the online continuous-flow enzymic reaction, using a fluorescent probe. The probe consists of fluorophore-labeled phosphate-binding protein, which shows a strong fluorescence enhancement upon binding to inorg. phosphate. To detect very small changes of the phosphate concn. in a postcolumn enzymic reaction medium, the enzymic removal of phosphate impurities from solvents, reagents, and samples was optimized for application in continuous flow. The potential of the

phosphate probe is demonstrated by monitoring the enzymic activity, i.e., the phosphate release, from alk. phosphatase. The selectivity of the phosphate readout, necessary to distinguish between phosphate contg. substrate or product and free inorg. phosphate released after enzymic conversion, is shown. The applicability of LC coupled to the enzymic assay using the phosphate readout was demonstrated by detection of tetramisole in a plant ext. as inhibitor of alk. phosphatase.

Parallel mass spectrometry allowed the simultaneous confirmation of the identity of the inhibitor.

RE.CNT 19 THERE ARE 19 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L10 ANSWER 2 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 2002:949675 CAPLUS
DN 139:134016
TI TEM Imaging of Mass-selected Polymer Molecules
AU Nasibulin, Albert G.; Kauppinen, Esko I.; Thomson, Bruce A.; Fernandez de la Mora, J.
CS Aerosol Technology Group, VTT Processes, FIN-02044, Finland
SO Journal of Nanoparticle Research (2002), 4(5), 449-453
CODEN: JNARFA; ISSN: 1388-0764
PB Kluwer Academic Publishers
DT Journal
LA English
AB Polyethylene glycol (PEG) mols. with masses below 1300 amu are electrosprayed (ES) from soln., mobility-selected at high resoln. in a differential mobility analyzer (DMA), collected on a grid and imaged by transmission electron microscopy (ES-DMA-TEM). The DMA resolves individual n-mers, and selects only one out of the many present in the original sample. Ion identity is established from **parallel mass spectra** (ES-MS). The images reveal spherical particles 1.46 nm in diam., in good agreement with the known ion mass and bulk d. The DMA-selection technique opens new paths for the study of very small particles.

RE.CNT 21 THERE ARE 21 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L10 ANSWER 3 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 2002:688793 CAPLUS
DN 137:358560
TI Thermal Stability of Self-Assembled Monolayers: Influence of Lateral Hydrogen Bonding
AU Valiokas, Ramunas; Oestblom, Mattias; Svedhem, Sofia; Svensson, Stefan C. T.; Liedberg, Bo
CS Division of Applied Physics, Division of Chemistry, Department of Physics and Measurement Technology, Linköping University, Linköping, S-581 83, Swed.
SO Journal of Physical Chemistry B (2002), 106(40), 10401-10409
CODEN: JPCBPK; ISSN: 1520-6106
PB American Chemical Society
DT Journal
LA English
AB Temp.-programmed desorption (TPD) of self-assembled monolayers (SAMs) on Au is studied by using **parallel mass spectrometry** (MS) and IR reflection-absorption spectroscopy (IRAS). Monolayers formed by HS(CH₂)_n-OH (n = 18, 22) and HS(CH₂)₁₅-CONH-(CH₂CH₂O)-H (EG1) are compared to reveal the influence of specifically introduced hydrogen-bonding groups on their thermal stability. The overall desorption process of the above mols. occurs in 2 main steps; a disordering of the alkyl chains followed by a complex series of decompn./desorption reactions. The final step of the process involves desorption of S from different chemisorption states. The amide-group-contg. SAM, which is stabilized by lateral hydrogen bonds, displays a substantial delay of the alkyl chain disordering by .apprx.50

K, as compared to the linear chain alcs. HS(CH₂)_n-OH. also, the decompn. of the alkyls and the onset of S desorption occur at a temp. that is higher by .apprx.25 K as compared to the HS(CH₂)₁₈-OH SAM. The desorption process is also studied for 2 oligo(ethylene glycol)-terminated SAMs, HS(CH₂)₁₅-X-(CH₂CH₂O)₄-H (EG4-SAMs), where X is -CONH- and -COO- linking groups. In addn. to the mol. chain disordering, the decompn./desorption process of the EG4-SAMs occurs in 2 steps. The 1st is assocd. with the loss of the oligomer portion and the 2nd with the desorption of the alkythiolate part of the mol. Study points out that lateral hydrogen bonding, introduced via amide groups, is a convenient way to improve the thermal stability of alkanthiolate SAMs.

RE.CNT 52 THERE ARE 52 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L10 ANSWER 4 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN

AN 2002:59880 CAPLUS

DN 136:262794

TI Neutral products from gas phase rearrangements of simple carbocations

AU Morton, Thomas Hellman

CS Department of Chemistry, University of California, Riverside, CA,
 92521-0403, USA

SO Advances in Gas Phase Ion Chemistry (2001), 4, 213-256

 CODEN: AGPCER; ISSN: 1071-9687

PB JAI Press Inc.

DT Journal; General Review

LA English

AB A review; analyzing the neutral products from ionic reactions in the gas phase provides information that cannot be gained by mass spectrometric methods alone. Neutrals have been recovered using three general techniques for generating ions in sufficient quantities: nuclear decay of multiply tritiated precursors, .gamma.-radiolysis studies, and electron bombardment flow (BBFlow) expts. Analyses of the uncharged reaction products of ion-mol. reactions are most effectively interpreted in conjunction with **parallel mass spectrometric** investigations. Taken together, these combined studies demonstrate the propensity of gaseous cations to undergo similar sorts of isomerizations as have been reported in condensed media. The absence of solvent and counterions makes it possible to produce ions in the gas phase that cannot be formed in soln. Despite the difference in reaction medium, the same two general categories of rearrangement-ring closure/ring opening and atom/group transfer-account for the variety of ion structures that give rise to the obsd. neutral products.

RE.CNT 95 THERE ARE 95 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L10 ANSWER 5 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN

AN 2000:106149 CAPLUS

DN 133:101503

TI A multiple electrospray interface for **parallel mass spectrometric** analyses of compound libraries

AU Wang, T.; Zeng, L.; Cohen, J.; Kassel, Daniel B.

CS CombiChem, Inc., San Diego, CA, USA

SO Combinatorial Chemistry and High Throughput Screening (1999), 2(6),
 327-334

 CODEN: CCHSFU; ISSN: 1386-2073

PB Bentham Science Publishers

DT Journal

LA English

AB A parallel spray interface for mass spectrometry is described. This new electrospray interface enables effluent flow streams from an array of HPLC columns to be sampled independently and sequentially on a chromatog. time-scale. Unlike our previously reported parallel LC-MS interface, which incorporated a dual-sheath spray interface accommodating up to four flow streams that are sampled simultaneously, this new interface permits

up to four columns to be sampled sequentially by means of a stepping motor and rotating plate assembly. Effluent flow streams from an array of four HPLC columns are connected to a parallel arrangement of electrospray needles co-axial to the mass spectrometer entrance aperture. Within the needle assembly, the individual spray tips are oriented in a circular array, where each needle is situated 90 degrees relative to one another for four-column operation. An eight-column system is described with needles spaced at 45 degree intervals. In between the needle assembly and the mass spectrometer entrance aperture is a Teflon disk with a through-hole that is mounted to a stepping motor assembly. By precisely controlling the stepping of the motor assembly, it is possible to "sample" each of the spray positions multiple times per s. By operating the quadrupole mass spectrometer in the single ion monitoring (SIM) mode, it was possible to acquire data at each of the spray positions during the course of the elution of compds. from the HPLC column array while maintaining both good sensitivity and peak shape. Preliminary results suggest this technique will be useful for high throughput combinatorial library anal. and profiling.

RE.CNT 17 THERE ARE 17 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L10 ANSWER 6 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN

AN 1998:148610 CAPLUS

DN 128:267854

TI Dual **parallel mass spectrometers** for
analysis of sphingolipid, glycerophospholipid and plasmalogen molecular
species

AU Byrdwell, Wm. Craig

CS PQS, NCAUR, ARS, USDA, Peoria, IL, 61604, USA

SO Rapid Communications in Mass Spectrometry (1998), 12(5), 256-272

CODEN: RCMSEF; ISSN: 0951-4198

PB John Wiley & Sons Ltd.

DT Journal

LA English

AB Anal. of phospholipids was performed using a liq. chromatog. sepn. with two mass spectrometers in parallel providing electrospray ionization (ESI) and atm. pressure chem. ionization (APCI) data simultaneously from a triple quadrupole instrument and a single quadrupole instrument, resp. The output from UV-Vis and evaporative light scattering detectors were also acquired by the two mass spectrometers, resp., for four detectors overall. This arrangement was used to identify and calc. area percents for mol. species of dihydrosphingomyelin (DHS) and sphingomyelin (SPM) in com. available bovine brain SPM, in human plasma ext. and in porcine lens ext. Mol. species of phosphatidylethanolamine and its plasmalogen, and phosphatidylcholine and its plasmalogen were identified and semi-quant. anal. performed. Com. available bovine brain SPM was found to contain 11.5% DHS and 88.5% SPM. The only DHS mol. species identified in human plasma was 16:0-DHS, at or below 1% of the sphingolipid content. Porcine lens membranes were found to contain 14.4% DHS and 85.6% SPM. Other findings reported here include: (1) phospholipids were found to undergo dimerization in the electrospray source, giving masses representing combinations of species present. (2) Triacylglycerols gave usable mass spectra under electrospray ionization conditions, as well as under APCI-MS conditions. (3) Triacylglycerols gave ammonium adducts as base peaks in their APCI mass spectra, which reduced fragmentation and increased the proportions of mol. ions. (4) Mass spectra were obtained for phospholipids which underwent both protonation and sodium adduct formation in different chromatog. runs. This paper was prepd. under the auspices of the US Government and it is therefore not subject to copyright in the US.

RE.CNT 26 THERE ARE 26 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L10 ANSWER 7 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN

AN 1998:45137 CAPLUS

DN 128:75016
TI Methanol Oxidation on Rhodium As Probed by Surface-Enhanced Raman and Mass Spectroscopies: Adsorbate Stability, Reactivity, and Catalytic Relevance
AU Williams, Christopher T.; Takoudis, Christos G.; Weaver, Michael J.
CS School of Chemical Engineering and Department of Chemistry, Purdue University, West Lafayette, IN, 47907, USA
SO Journal of Physical Chemistry B (1998), 102(2), 406-416
CODEN: JPCBFK; ISSN: 1089-5647
PB American Chemical Society
DT Journal
LA English
AB

The relationship between surface speciation and catalytic activity/selectivity during methanol oxidn. on polycryst. rhodium under ambient-pressure flow-reactor conditions was studied from 25 to 500 .degree.C by means of surface-enhanced Raman spectroscopy (SERS) along with parallel mass spectrometric (MS) measurements. By utilizing SERS-active Rh films formed by electrodeposition onto gold, the former technique provides in situ surface vibrational spectra with unique sensitivity under these demanding conditions, enabling adsorbed species to be probed in real time (.apprxeq.1 s) for comparison with the overall kinetics as evaluated by MS. Exposure of Rh to O2-free methanol yielded no detectable vibrational bands between 25 and 500.degree., although methanol decompn. to form CO and H2 was evident from MS. The presence of even subunity molar ratios of oxygen, however, yielded rich SER spectra, highlighted by bands indicative of CO(ads) (.nu.Rh-CO = 465 cm-1, .nu.Rh-CO = 2000 cm-1). The catalytic selectivity toward CO2 (vs. CO) gaseous product formation decreased markedly around the desorption temp. of CO(ads) .apprxeq. 350.degree. under these conditions. This is consistent with the facilitation of CO2 prodn. by the presence of CO(ads). Complete selectivity toward exhaustive methanol oxidn. (i.e., CO2, H2O formation) was obsd. in oxygen-rich methanol mixts., adsorbed CO now being absent at all temps. The CO2 prodn. occurs partly via methanolic C-O cleavage as deduced by 18O2 substitution. The presence of rhodium oxide (Rh2O3) was diagnosed with such reactant mixts. above ca. 300 .degree.C from the characteristic 500-580 cm-1 .nu.Rh-O bands. The kinetics of formation and removal of the oxide were probed by gas flow switching coupled with transient SERS measurements. The oxide formation rates following O2 exposure are slowed markedly (>100-fold) by the presence of even a small (5%) methanol mole fraction. Switching to pure methanol results in very rapid oxide redn., so that, for example, removal is complete within ca. 1s at 350.degree. with 100 Torr of CH3OH. Examn. of the transient oxide removal kinetics as a function of temp. and methanol pressure revealed a transition from strongly activated to essentially T-independent behavior at lower pressures and/or higher temps. This is indicative of a change from rate-detcg. removal of oxygen from the oxide lattice to a subsequent step involving formation of and/or reaction with an adsorbed methanol scavenger. While such reactivity earmarks the oxide as a potential reaction intermediate, the overall catalytic turnover rates for methanol oxidn. are nonetheless faster than can be accommodated on this basis.

RE.CNT 45 THERE ARE 45 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L10 ANSWER 8 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 1997:410724 CAPLUS
DN 127:155996
TI Probing Combinatorial Library Diversity by Mass Spectrometry
AU Demirev, Plamen A.; Zubarev, Roman A.
CS Division of Ion Physics The Aangstroem Laboratory, Uppsala University, Uppsala, S-751 21, Swed.
SO Analytical Chemistry (1997), 69(15), 2893-2900
CODEN: ANCHAM; ISSN: 0003-2700
PB American Chemical Society
DT Journal

LA English
AB The feasibility of a massively **parallel mass spectrometric** method for probing combinatorial library diversity is addressed theor. for the example of computer-generated mass distributions of combinatorially synthesized peptide libraries contg. between two and seven amino acids. The authors study the behavior of several global (integral) parameters of such mass distributions-mass centroid, dispersion, skewness, and kurtosis. The centroid and dispersion carry information that may characterize the completeness of the synthetic effort. Local mass distribution parameters, e.g., mass d. (no. of peptides per mass interval), are also examd. The practical implementation and eventual limitations of such an approach are discussed as well.

L10 ANSWER 9 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN

AN 1996:667782 CAPLUS

DN 125:320017

TI Detection of the picolinic acid biomarker in *Bacillus* spores using a potentially field-portable pyrolysis-gas chromatography-ion mobility spectrometry system

AU Snyder, A. Peter; Thornton, Sidney N.; Dworzanski, Jacek P.; Meuzelaar, Henk L. C.

CS Dev. Eng. Cent., U.S. Army Edgewood Res., Aberdeen Proving Ground, MD, 21010-5423, USA

SO Field Analytical Chemistry and Technology (1996), 1(1), 49-59

CODEN: FACTFR; ISSN: 1086-900X

PB Wiley

DT Journal

LA English

AB The absence of a field-portable device to provide real-time detection of Gram-pos. bacterial spores has prompted the interfacing of a pyrolysis (Py) module to an existing, hand-held gas-chromatog.-ion-mobility spectrometry (GC/IMS) device. In this configuration, spore detection is achieved by the observation of picolinic acid (PA), which is the most characteristic pyrolysis decompn. product of the parent dipicolinic (2,6-pyridinedicarboxylic) acid (DPA). Pos. identification of PA was demonstrated using a lab.-based GC instrument with dual, **parallel mass spectrometry** (MS) and IMS detectors. Spores and vegetative microorganisms of the genus *Bacillus* were characterized by the presence and absence of DPA, resp., and the picolinic acid marker was identified from the GC/IMS and GC/MS profiles. A field-portable prototype Py-GC/IMS system is described and appears to provide similar bioanal. information with respect to the lab.-based system. Preliminary results of this study indicate that the degree of compd. sepn. afforded by a short GC capillary column guards against common environmental interferences including urban particulate matter and biol. particles such as fungal spores and pollen.

L10 ANSWER 10 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN

AN 1987:429781 CAPLUS

DN 107:29781

TI Determination of the nuclear reactor burning process balance by gamma spectrometry of fission products. Part V - Determination of the isotopic composition of irradiated uranium

AU Bulovic, V.; Maksimovic, Z.; Krtil, J.; Sus, F.; Kloseva, E.

CS Boris Kidric Inst. Nucl. Sci., Vinca, Yugoslavia

SO Jaderna Energije (1987), 33(1), 8-11

CODEN: JADEAQ; ISSN: 0448-116X

DT Journal

LA English

AB The possibility of detg. the isotopic compn. of irradiated U fuel of a heavy water reactor on the basis of .gamma.-spectrometry of fission products was exptl. tested. The testing was performed upon spent fuel from unenriched U. For detg. the fission products (¹⁰⁶Ru, ¹³⁴Cs and ¹³⁷Cs) a spectrometer with a Ge(Li) detector was used. The accuracy of

the results obtained for the compn. of U was tested through its
parallel mass-spectrometric analyses.

L10 ANSWER 11 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 1987:112632 CAPLUS
DN 106:112632
TI The thermal decomposition of strontium fluorophosphate hydrate
(SrPO3F.cntdot.H2O)
AU Menz, D. H.; Heide, K.; Kunert, C.; Mensing, C.; Kolditz, L.
CS Zentralinst. Anorg. Chem., Dtsch. Akad. Wiss., Berlin, DDR-1199, Ger. Dem.
Rep.
SO Zeitschrift fuer Anorganische und Allgemeine Chemie (1986), 540-541, 191-7
CODEN: ZAACAB; ISSN: 0044-2313
DT Journal
LA German
AB The thermal decompn. of SrPO3F.H2O was studied by complex thermal anal.
The thermogravimetric study was completed by simultaneous and
parallel mass spectrometric anal. of the gas
phase. During the 1st state of thermal decompn. .apprx.0.8 mol water is
lost. Then a partial hydrolysis takes place and HF is formed. The
formation of POF3 is a multistage mechanism without effect of H2O at
>500.degree.. The partial reactions leading to .alpha.-Sr2P2O7 and SrF2
>600.degree. and to .alpha.-Sr2P2O7 and Sr5(PO4)3F >750.degree. were
formulated and the exptl. and calcd. mass loss were compared.

L10 ANSWER 12 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 1985:196916 CAPLUS
DN 102:196916
TI Thermal decomposition of calcium phosphorofluoridate dihydrate
(CaPO3F.2H2O)
AU Heide, K.; Menz, D. H.; Schmidt, C.; Kolditz, L.
CS Sekt. Chem., Friedrich-Schiller-Univ., Jena, DDR-6900, Ger. Dem. Rep.
SO Zeitschrift fuer Anorganische und Allgemeine Chemie (1985), 520, 32-8
CODEN: ZAACAB; ISSN: 0044-2313
DT Journal
LA German
AB The thermal decompn. of CaPO3F.2H2O was studied by thermogravimetry under
inert conditions. A **parallel mass
spectrometric** anal. of gases produced was made. With the use of
an effusion cell a quasiequil. evapn. in the vicinity of the ion source of
the spectrometer was achieved. The results are comparable with the
thermogravimetric anal. under normal pressure. During 1st stage of
thermal decompn. 1 mol H2O was lost. The further course is detd. by
release of HF and POF3. The several steps of decompn. leading to
.alpha.-Ca2P2O7 at >360.degree. are discussed.

L10 ANSWER 13 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 1985:178482 CAPLUS
DN 102:178482
TI Method and apparatus for **parallel mass
spectrometry**
PA Chang, Chuang, USA
SO Jpn. Kokai Tokkyo Koho, 16 pp.
CODEN: JKXXAF
DT Patent
LA Japanese
FAN. CNT 1

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 59176663	A2	19841006	JP 1984-41208	19840303
US 4507555	A	19850326	US 1983-472161	19830304
PRAI US 1983-472161		19830304		

AB The design is claimed of a **parallel mass
spectrometric** app. joined in tandem with a gas chromatograph.

L10 ANSWER 14 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 1985:159764 CAPLUS
DN 102:159764
TI **Parallel mass spectrometry** for high
performance GC and LC detection
AU Chang, C.
CS Wright State Univ., Dayton, OH, USA
SO American Laboratory (Shelton, CT, United States) (1985), 17(3), 59-64, 66
CODEN: ALBYBL; ISSN: 0044-7749
DT Journal; General Review
LA English
AB A review with 15 refs. Problem areas in using conventional scanning mass spectrometers for high-performance gas chromatog. (GC) and liq. chromatog. (LC) detection are discussed. The potential use of **parallel mass spectrometers** to avoid these problems is also discussed.

L10 ANSWER 15 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 1976:600351 CAPLUS
DN 85:200351
TI Valence level photoelectron spectra of some heavy group 4-6 diatomic molecules
AU Wu, M.; Fehlner, T. P.
CS Dep. Chem., Univ. Notre Dame, Notre Dame, IN, USA
SO Journal of the American Chemical Society (1976), 98(24), 7578-85
CODEN: JACSAT; ISSN: 0002-7863
DT Journal
LA English
AB The He I photoelectron spectra of GeS, GeSe, SnS, SnTe, and PbTe in the gas phase were obtained by the photoionization of the vapors above appropriate solids at 700-1000.degree.K. Spectra are assigned by using obsd. relative band areas, vibrational fine structure, and spin-orbit splitting along with electron impact ionization potentials and **parallel mass spectrometric** studies. There is significant mixing of the .SIGMA.1/2 and .PI.1/2 states in the heavier species. Distinct differences between the .PI. states of light and heavy diatomics are obsd. Similarities and differences between the valence regions of group 4-6 diatomics and diatomics of group 5-5 and group 3-7 are also reported.

L10 ANSWER 16 OF 16 CAPLUS COPYRIGHT 2003 ACS on STN
AN 1959:15979 CAPLUS
DN 53:15979
OREF 53:2923h-i,2924a
TI Anomalous behavior of gem-diethers in the mass spectrometer
AU LeBlanc, R. Bruce
CS Dow Chem. Co., Freeport, TX
SO Anal. Chem. (1958), 30, 1797-9
CODEN: ANCHAM; ISSN: 0003-2700
DT Journal
LA Unavailable
AB gem-Diether (compds. with 2 alkoxy groups on the same C atom) were measured in the mass spectrometer. They give different spectra, depending on whether the filament is bare W or carbonized W. The carbonized filament gives the normal spectrum. The bare filament causes a partial decompn. into a vinyl ether and alc. For example, MeCH(OEt)2 on decompn. yields EtOH and CH2:CHOEt. For reliable analysis of the compds. a carbonized filament is recommended.

(FILE 'WPIDS' ENTERED AT 16:18:38 ON 24 DEC 2003)

FILE 'USPATFULL' ENTERED AT 16:21:50 ON 24 DEC 2003

L4 1 S (PARALLEL(W) MASS (W)SPECTROMETRY)/TI,AB,CLM
L5 2 S (PARALLEL(W) MASS (W)SPECTROMETRY)/TI,AB,CLM
L6 3 S (PARALLEL(W) MASS (W)SPECTR?)/CLM,AB,TI
L7 3 S (PARAL?(W) MASS (W)SPECTR?)/CLM,AB,TI
L8 3 S (PARAL?(W) MASS (W)SPECTR?)
L9 8 S L8 AND PROTEIN

=> d bib,kwic 1-8

L9 ANSWER 1 OF 8 USPATFULL on STN
AN 2003:173349 USPATFULL
TI System and method for high throughput screening of droplets
IN Hess, Robert, Arlington, MA, UNITED STATES
Brenan, Colin, Marblehead, MA, UNITED STATES
Linton, John, Lincoln, MA, UNITED STATES
Ozbal, Can, Cambridge, MA, UNITED STATES
Green, Donald, Watertown, MA, UNITED STATES
Hunter, Ian, Lincoln, MA, UNITED STATES
PI US 2003119193 A1 20030626
AI US 2002-267912 A1 20021008 (10)
RLI Continuation-in-part of Ser. No. US 2001-842361, filed on 25 Apr 2001,
PENDING
DT Utility
FS APPLICATION
LREP BROMBERG & SUNSTEIN LLP, 125 SUMMER STREET, BOSTON, MA, 02110-1618
CLMN Number of Claims: 90
ECL Exemplary Claim: 1
DRWN 14 Drawing Page(s)
LN.CNT 2283
CAS INDEXING IS AVAILABLE FOR THIS PATENT.
SUMM . . . speeds at which large numbers of samples can be analyzed.
Unlike optical-based assays in which samples can be analyzed in
parallel, mass spectrometry is a serial
process in which sample must be analyzed one-at-a-time. Typically, a
slow desalting step or purification step is. . .
SUMM . . . seconds. In various embodiments, the rate is substantially one
assay per second. The reaction may also be buffered only by
proteins intrinsic to the assay such as the enzyme in an enzyme
inhibition assay.
DETD [0121] .alpha.-Chymotrypsin is a protease that cleaves **proteins**
and peptides at aromatic amino acids such as phenylalanine, tyrosine,
and tryptophan. The example assay attempts to discover inhibitors of. . .

L9 ANSWER 2 OF 8 USPATFULL on STN
AN 2002:268969 USPATFULL
TI Mass spectrometer apparatus for analyzing multiple fluid samples
concurrently
IN Moini, Mehdi, Austin, TX, United States
Jiang, Longfei, Austin, TX, United States
PA Board of Regents, The University of Texas System, Austin, TX, United
States (U.S. corporation)
PI US 6465776 B1 20021015
AI US 2000-586588 20000602 (9)
DT Utility
FS GRANTED
EXNAM Primary Examiner: Lee, John R.; Assistant Examiner: Vanore, David A.
LREP Fulbright & Jaworski, LLP
CLMN Number of Claims: 24
ECL Exemplary Claim: 1

DRWN 2 Drawing Figure(s); 2 Drawing Page(s)

LN.CNT 713

CAS INDEXING IS AVAILABLE FOR THIS PATENT.

SUMM . . . of dual ESI sprayers have been tried with a Y-shaped orifice defined within the nozzle in order to investigate electrosprayed **proteins** using ion-ion or ion-molecule reactions. In particular the accurate measurement of masses of organic compounds has been another use of. . .

SUMM . . . Corporation on a "Hot Gas Sampling"; and U.S. Pat. No. 4,507,555 patented Mar. 26, 1985 to C. Chang on a "**Parallel Mass Spectrometer**"; and U.S. Pat. No. 4,562,351 patented Dec. 31, 1985 to P. Atherton et al and assigned to VG Instruments Group. . .

L9 ANSWER 3 OF 8 USPATFULL on STN

AN 2002:191521 USPATFULL

TI Massive parallel method for decoding DNA and RNA

IN Ju, Jingyue, Englewood Cliffs, NJ, UNITED STATES

Li, Zengmin, New York, NY, UNITED STATES

Edwards, John Robert, New York, NY, UNITED STATES

Itagaki, Yasuhiro, New York, NY, UNITED STATES

PI US 2002102586 A1 20020801

US 6664079 B2 20031216

AI US 2001-972364 A1 20011005 (9)

RLI Continuation-in-part of Ser. No. US 2000-684670, filed on 6 Oct 2000, PENDING

PRAI US 2001-300894P 20010626 (60)

DT Utility

FS APPLICATION

LREP John P. White, Cooper & Dunham LLP, 1185 Avenue of the Americas, New York, NY, 10036

CLMN Number of Claims: 60

ECL Exemplary Claim: 1

DRWN 28 Drawing Page(s)

LN.CNT 1869

CAS INDEXING IS AVAILABLE FOR THIS PATENT.

SUMM [0034] The invention provides a **parallel mass spectrometry** system, which comprises a plurality of atmospheric pressure chemical ionization mass spectrometers for parallel analysis of a plurality of samples. . .

DRWD [0057] FIG. 24: **Parallel mass spectrometry** system for DNA sequencing. Example shows three mass spectrometers in parallel. Samples are injected into the ion source where they. . .

DETD . . . In one embodiment, the mass tag is a 2-nitro-.alpha.-methyl-3,4-dimethoxybenzyl group. In one embodiment, the mass tag is detected using a **parallel mass spectrometry** system which comprises a plurality of atmospheric pressure chemical ionization mass spectrometers for parallel analysis of a plurality of samples. . .

DETD [0133] The invention provides a **parallel mass spectrometry** system, which comprises a plurality of atmospheric pressure chemical ionization mass spectrometers for parallel analysis of a plurality of samples. . .

DETD . . . areas of biomedical research. Though these ionization methods are suitable for the analysis of bioorganic molecules, such as peptides and **proteins**, improvements in both detection and sample preparation are required for implementation of mass spectrometry for DNA sequencing applications. Since the. . .

DETD [0151] The photocleavable 2-nitrobenzyl moiety has been used to link biotin to DNA and **protein** for efficient removal by UV light (.about.350 nm) (Olejnik et al. 1995, 1999). In the approach disclosed herein the 2-nitrobenzyl. . .

DETD . . . not capped. As discussed above, the photo cleavable 2-nitro benzyl moiety has been used to link biotin to DNA and **protein** for efficient removal by UV light (.about.350 nm) irradiation (Olejnik

et al. 1995, 1999). Four different 2-nitro benzyl groups with. . . .
 DETD [0172] To make mass spectrometry competitive with a 96 capillary array
 method for analyzing DNA, a **parallel mass**
spectrometer approach is needed. Such a complete system has not
 been reported mainly due to the fact that most of the. . . .
 DETD [0173] A complete **parallel mass spectrometry**
 system includes multiple APCI sources interfaced with multiple
 analyzers, coupled with appropriate electronics and power supply
 configuration. A mass spectrometry. . . figures show a system with
 three mass spectrometers in parallel. Higher throughput is obtained
 using a greater number of in **parallel mass**
spectrometers.
 CLM What is claimed is:
 26. The method of claim 17, wherein the mass tag is detected using a
parallel mass spectrometry system which
 comprises a plurality of atmospheric pressure chemical ionization mass
 spectrometers for parallel analysis of a plurality of samples. . . .
 54. A **parallel mass spectrometry** system,
 which comprises a plurality of atmospheric pressure chemical ionization
 mass spectrometers for parallel analysis of a plurality of samples. . .

L9 ANSWER 4 OF 8 USPATFULL on STN
 AN 2002:133513 USPATFULL
 TI Proteomic analysis by **parallel mass**
spectrometry
 IN LaDine, James R., Uxbridge, MA, UNITED STATES *ap 5*
Jardine, Ian, Los Gatos, CA, UNITED STATES
Storv, Mike S., Los Gatos, CA, UNITED STATES
 PI US 2002068366 A1 20020606
 AI US 2001-835273 A1 20010413 (9)
 PRAI US 2000-196889P 20000413 (60)
 DT Utility
 FS APPLICATION
 LREP John J. Gagel, Fish & Richardson P.C., 225 Franklin Street, Boston, MA,
 02110-2804
 CLMN Number of Claims: 28
 ECL Exemplary Claim: 1
 DRWN 7 Drawing Page(s)
 LN.CNT 1181
 TI Proteomic analysis by **parallel mass**
spectrometry
 SUMM [0002] This invention relates to proteomic analysis by **parallel**
mass spectrometry.
 SUMM [0003] Within a typical cell there are several thousand **proteins**
 , its "proteome," which carry out the metabolic work of the cell. These
proteins are in constant interplay with one another, and with
 every other sort of biomolecule found within a cell. The
proteins physically interact, or bind, to each other and to
 common secondary molecules. The result of such interactions is a fine
 control and balancing of metabolic functions. For example, one
protein may increase or decrease the function of another
protein by binding to it and altering its structure by the
 addition or removal of a modifying group such as a phosphate. Another
 mode of action is for one **protein** to produce more or less of a
 secondary substance that interacts allosterically with a second
protein (or multiple second **proteins**) to modulate its
 function. Analysis of the abundance of **proteins** can therefore
 be useful in elucidating the molecular basis of differences brought
 about by diseases or by therapeutic treatments
 SUMM [0004] A number of techniques have been suggested for analyzing cellular
proteins, including, for example, two-dimensional
 electrophoresis followed by mass spectrometry. In the case of
 two-dimensional electrophoresis, a **protein** sample is placed in

array to a common computing device, said mass spectral data being indicative of the identity and the abundance of **protein** in said multiple sample, and correlating said mass spectral data as a function of time.

2. The method of claim 1 comprising displaying said correlated data as a function of **protein** identity, **protein** abundance, and time.

4. The method of claim 1 comprising identifying **proteins** based on changes in abundance as a function of time.

6. The method of claim 4 comprising analyzing 500 **proteins** or more.

7. The method of claim 6 comprising analyzing 5000 **proteins** or more.

22. A method for analysis of **proteins** in a biological system comprising: providing a biological system containing **proteins**; exposing the biological system to a stimulus; after exposing the biological system to the stimulus, sampling the biological system at multiple time intervals to obtain multiple samples; treating the multiple samples by a separation technique to provide multiple **protein** samples suitable for analysis by mass spectrometry; providing a parallel array of mass spectrometer systems capable of simultaneous analysis of as many **protein** samples as there are spectrometer systems in said array; analyzing the multiple **protein** samples in said parallel array of mass spectrometry systems to generate mass spectral data indicative of the identity and the abundance of **proteins** in said multiple **protein** samples; and in a common electronic computing device communicating with each of said mass spectrometry systems, correlating said mass spectral.

27. The system of claim 26 wherein the analysis includes analysis of about 500 **proteins** or more.

L9 ANSWER 5 OF 8 USPATFULL on STN

AN 1999:18911 USPATFULL

TI Methods and apparatus for sequencing polymers with a statistical certainty using mass spectrometry

IN Patterson, Dale H., Nashua, NH, United States

PA PerSeptive Biosystems, Inc., Framingham, MA, United States (U.S. corporation)

PI US 5869240 19990209

AI US 1995-447175 19950519 (8)

DT Utility

FS Granted

EXNAM Primary Examiner: Jones, W. Gary; Assistant Examiner: Rees, Dianne

LREP Testa, Hurwitz & Thibault, LLP

CLMN Number of Claims: 47

ECL Exemplary Claim: 1

DRWN 16 Drawing Figure(s); 16 Drawing Page(s)

LN.CNT 1668

CAS INDEXING IS AVAILABLE FOR THIS PATENT.

SUMM . . . complete primary structure identification. To date, Edman sequencing and adaptations thereof are the most widely used tools for sequencing certain **protein** and peptides residue by residue, while the enzymatic synthesis method is preferred for sequencing oligonucleotides.

SUMM In the case of **protein** and peptide sequencing, C-terminal sequencing via chemical methods has proven particularly difficult while being only marginally effective, at best. (See, e.g., Spiess, J. (1986)

Methods of **Protein** Characterization: A Practical Handbook (Shively, J. E. ed., Humana Press, N.J.) pp. 363-377; Tsugita et al. (1994) J. **Protein** Chemistry 13:476-479). Consequently, the C-terminus remains a region often not analyzed because of lack of a dependable method.

SUMM offer a simple approach by which amino acids can be sequentially cleaved residue by residue from the C-terminus of a **protein** or a peptide. Carboxypeptidase Y (CPY), in particular, is an attractive enzyme because it non-specifically cleaves all residues from the

SUMM by residue. Not only is this approach labor-intensive, but it is complicated by amino acid contaminants in the enzyme and **protein**/peptide solutions, as well as by enzyme autolysis. A further hindrance to any sequencing effort of this type is the absolute.

SUMM analysis such as field desorption (Hong et al. (1983) Biomed. Mass Spectrom. 10:450-457), electrospray (Smith et al. (1993) 4 Techniques **Protein** Chem. 463-470), and thermospray (Stachowiak et al. (1988) J. Am. Chem. Soc. 110:1758-1765), it is possible to perform direct mass. . . .

SUMM digestion of peptides has been combined with other mass spectrometry methods such as plasma desorption (Wang et al. (1992) Techniques **Protein** Chemistry III (ed., R. H. Angeletti; Academic Press, N.Y.) pp. 503-515).

SUMM obtaining sequence information that incorporates a data interpretation strategy based on integrating mass-to-charge ratio data obtained from a plurality of **parallel mass spectra**.

SUMM The claimed methods are applicable to any polymer, including biopolymers such as DNAs, RNAs, PNAs, **proteins**, peptides and carbohydrates, and modified forms of these polymers. The set of polymer fragments may be created by hydrolysis of. . . .

DETD moiety. In a currently preferred embodiment, the polymer is a biopolymer selected from, but not limited to, the following group: **proteins**, peptides, DNAs, RNAs, PNAs (peptide nucleic acids), carbohydrates and modified forms thereof.

DETD The claimed invention can be applied to the sequencing of any natural biopolymer such as **proteins**, peptides, nucleic acids, carbohydrates, etc., as well as synthetic biopolymers such as PNA and phosphotiolated nucleic acids. The ladders could. . . .

DETD collective truncated hydrolyzed polymer fragments. In this manner, for example, sequence information relating to the amino acid sequence of a **protein** can be obtained using carboxypeptidase Y, an agent which acts at the carboxy terminus. By using the methods disclosed herein to generate a series of **protein** hydrolysates related one to the other by consecutive, repetitive liberation of amino acid residues, the skilled artisan can reconstruct the primary sequence of the intact **protein** polymer as described in further detail below.

DETD invention for this purpose. Thus the above-described subtractive-type sequencing method, through which repetitive removal of successive amino-terminal residues from a **protein** polymer can occur, can also be accomplished with hydrolyzing agents other than enzymes as disclosed herein.

DETD As disclosed herein, this strategy can be applied to the sequencing of any natural biopolymer such as **proteins**, peptides, nucleic acids, carbohydrates, etc. as well as synthetic biopolymers such as PNA and phosphotiolated nucleic acids. The ladders can. . . .

CLM What is claimed is:

5. The method of claim 4 wherein the biopolymer is selected from the group consisting of DNAs, RNAs, PNAs, **proteins**, peptides, carbohydrates and modified forms thereof.

22. The method of claim 21 wherein the biopolymer is selected from the